

level is kept below the level of the spout. To lower the water level in a SOCK, the diaphragm can be pumped with the index finger and jets of water will be ejected from the spout and no specimens, if they are crawling mollusks. This apparatus was made specifically for collecting sacoglossan and nudibranch mollusks in the size ranges of < 2.5 cm. The O-ring and medium gauge dental latex can easily be substituted by a string and a piece of toy balloon or rubber glove. If you keep a few "Acadian SOCKs" in your car, then a complete collecting kit for sea shore or lake is always at hand and the expendable parts can be purchased or replaced afresh in the shopping center of even the smallest town.

My tide pool technique was to carry a basket loaded with 10 Acadian SOCKs. Each species of nudibranch collected was kept in a separate SOCK. If occasional nudibranchs were too large or some of the dorids too stiff, then the rubber was flipped off and the specimen dropped in with the others and the diaphragm replaced. Specimens to be collected must be submerged, but it does not matter if the SOCK has air or water or both in it. The SOCK can be partly in air and water or totally submerged. Specimens were scraped off the underside of overturned rocks with the SOCK spout (or forceps) and positioned in the aperture of the spout. The diaphragm was depressed and the spout submerged in the nearest pool of water, releasing the finger pressure at the same time. The mucoid blob from the rock surface, now submerged in sea water in the SOCK, would then assume its natural shape and was examined with a hand lens.

The principle employed in the Acadian SOCK apparatus is applicable to larger and smaller diameter spouts and cylinders. As long as the volume displaced by depressing the rubber lid is greater than that of the spout, then any small or delicate organism can be drawn back into the cylinder within a gentle stream of water. I have field-tested two sizes of SOCK apparatus in 7 countries since January 1968 and can attest to their effectiveness. In fact, I have had the pleasure of using SOCK apparatus with SCUBA gear in *Zostera* beds and in algal jungles. One need only shake specimens off the plants or epizoans and then leisurely "pick them out of the air" as they drift about; or pieces of plants or hydroids or polyzoans with attached nidosomes and adult nudibranchs can be broken off and sucked up together.

The uninterrupted hours necessary to evolve the design of the Acadian SOCK were made possible by a year of Sabbatical leave from Acadia University, and a Canadian National Research Council Senior Research Fellowship. The author is further indebted to Professor George Hughes, Chairman, Department of Zoology, Bristol University, England, for providing workshop facilities where the prototypes and present model were constructed.

Technique for Extraction and Mounting of Gastropod Radulae

BY

GEORGE E. RADWIN

Curator of Marine Invertebrates
San Diego Natural History Museum
San Diego, California 92112

I HAVE USED the following procedure for extracting and mounting radulae for a number of years. It is largely a pragmatic technique developed by a process of trial and error.

The extraction procedure differs depending on the size and condition of the specimen from which the radula is to be extracted. If a dried specimen is to be used, it should be boiled in a saturated solution of trisodium phosphate (TSP). If the animal has dried deeply retracted into the shell such boiling in TSP will usually cause softening and swelling sufficient to allow it to be seen and removed from the shell. If the specimen is small, the entire animal is usually treated. If a large animal is used, the first step must be to remove the buccal mass, lying immediately behind the foot and head areas and just below the mantle cavity. The softened animal (or the buccal mass) is then boiled in a concentrated solution of sodium hydroxide for a short time (3-10 seconds), reducing it to a viscous brown film on the surface of the NaOH solution. This brown scum is removed by pipette and transferred to a watch glass with 70% ethanol. Gently swirling the watch-glass, while carefully viewing it under a dissection microscope will generally cause the brown film to dissipate and leave the radula as a highly refractive filamentous object. The radula is generally a ribbon-like structure which may be removed by the aid of needles, and stored in a vial of alcohol. If the radula of a toxoglossate species [Conidae, Terebridae, Turridae] is extracted, the radula may be a packet of tiny "darts."

The radula under microscopic scrutiny is transferred to a depression slide containing a drop of eosin stain and allowed to remain there until the drop has dried. Using a disposable hypodermic syringe the depression slide is carefully filled with 70% ethanol. When the radula has been sufficiently destained (5-10 minutes), the 70% ethanol is removed, using the hypodermic syringe, and replaced with 95% ethanol. The lightly stained radula is allowed to remain in the 95% ethanol for about 5 minutes, then transferred with needles to a slide with a drop of xylene or water-based mounting medium. Still viewing the radula under the dissection microscope, a small section of the radula is removed intact and segregated, the remainder

of the radula being dissociated with a pair of very finely ground needles. This step precludes confusion due to overlapping teeth. If xylene is used, the radular material is covered with a drop or more of canada balsam mounting medium. If a water-based mounting medium is used, the step involving xylene must be omitted and the radular material is transferred directly from 95% ethanol to the mounting medium, and covered with a cover slip. The prepared slide may be viewed after 2 or 3 days, by which time the radular teeth have stopped rotating. For best results a compound microscope capable of a magnification of at least $400\times$ should be used.

MATERIALS NEEDED

1. saturated solution of trisodium phosphate
2. concentrated solution of sodium hydroxide
3. alcohol lamp, ringstand, ring and wire mesh
4. 3 disposable hypodermic syringes
5. 2 very finely ground wooden-handled dissecting needles
6. 70% ethanol
7. 95% ethanol
8. xylene
9. canada balsam or other permanent mounting medium
10. eosin or eosin-type stain
11. 3 depression slides
12. 1 watch glass (1 - 2 inches in diameter)
13. 2 microscopes (1 dissection, 1 compound)
14. microscope slides and coverslips

